

CRB Bore Hole Filter

December 7, 2016

1 Methods

The stoppers are 1.125" on larger side. Attached are 3 pictures - The charcoal filters are for filtering incoming air. Push these filters in so the thin part (shown with forceps) will go a bit into the hole, inside the chewed part of palm just behind the beetle. Then push insert white filters in on the opposite side (outside the plant and attach that with bungee cord or rope to the plant. Put the pump on the filter and connect to the battery. 24 h filtering should be enough. Once the aeration is done, replace rubber septum and little plug on each filter and put it back in the bag, wrap it well (heat seal the bag if possible) mark location, sex of the beetle, date. Store in freezer until all samples are done. You have 12 filters. 6 for each sex if possible.

nnn



Figure 1:



Figure 2:




Figure 3: